

	Document ID:	TDS-BBL-001-100ML	Version:	002
	Date of Issue:	10-JAN-2023	Approved by:	Dr. Iman Kamranfar
	Review Date:	10-JAN-2025	Signature:	
	Title:	TECHNICAL DATASHEET		

Sodium Bicarbonate Solution (7.5%)

Filtration, Treatment	0.1µm sterile filtered
Product Code	BBL-001-100ML
Shelf Life	4 Years from DOM
Storage Temperature	15 - 30 °C
Shipping Temperature	ambient

QC Specifications

Physical and Chemical Analysis	Method	Specifications	Units
Appearance	Visual	clear, colourless solution	n/a
pH at RT	Electronic pH Meter	7.8 - 8.5	n/a
Osmolality	Osmometer	Test and report	mOsm/kg
Endotoxin	LAL Kinetic	≤ 1.0	EU/ml
Sterility			
Aerobic Bacteria	Internally Validated	Not detected	n/a
Anaerobic Bacteria	Internally Validated	Not detected	n/a
Fungi (Yeast & Mold)	Internally Validated	Not detected	n/a
Mycoplasma	qPCR	Not detected	n/a

Formulation

Components	CAS number	Concentration (mg/L)
Sodium bicarbonate	144-55-8	75000.00

Product use: This product is intended for laboratory use only.