



Document ID:	TDS-MCL-045-500ML	Version:	004
Date of Issue:	30-APR-2021	Approved by:	Anna Becker
Review Date:	10-JAN-2022	Signature:	<i>Anna Becker</i>
Title:	<b>TECHNICAL DATA SHEET</b>		

<b>Product Name</b>	Dulbecco's Modified Eagle's Medium (DMEM)
<b>Filtration/Treatment</b>	With high Glucose, with stable Glutamine and 25 mM HEPES
<b>Product Codes</b>	MCL-045-500ML
<b>Shelf Life</b>	18 months
<b>Storage Temperature</b>	2 to 8 °C
<b>Shipping Temperature</b>	ambient
<b>Substance</b>	<b>Concentration (mg/L)</b>
<b>Amino Acids</b>	
Glycine	30.00
L-alanyl-L-glutamine	868.06
L-Arginine, Hydrochloride	84.00
L-Cystine, Dihydrochloride	63.00
L-Histidine, Hydrochloride, Monohydrate	42.00
L-Isoleucine	105.00
L-Leucine	105.00
L-Lysine, Hydrochloride	146.00
L-Methionine	30.00
L-Phenylalanine	66.00
L-Serine	42.00
L-Threonine	95.00
L-Tryptophan	16.00
L-Tyrosine, Disodium, Dihydrate	104.00
L-Valine	94.00
<b>Vitamins</b>	
Choline Chloride	4.00
D-Calcium Pantothenate	4.00
Folic Acid	4.00
Myo-Inositol	7.20
Niacinamide	4.00
Pyridoxine, Hydrochloride	4.00
Riboflavin	0.40
Thiamine, Hydrochloride	4.00
<b>Other Components</b>	
D-Glucose, Anhydrous	4500.00
HEPES	5958.00
Phenol Red, Sodium Salt	15.00
Pyruvic Acid, Sodium Salt	110.00



Document ID:	TDS-MCL-045-500ML	Version:	004
Date of Issue:	30-APR-2021	Approved by:	Anna Becker
Review Date:	10-JAN-2022	Signature:	
Title:	<b>TECHNICAL DATA SHEET</b>		

<b>Inorganic Salts</b>	<b>Concentration (mg/L)</b>
Calcium Chloride, Anhydrous	200.00
Ferric Nitrate, Nonahydrate	0.10
Magnesium Sulfate, Anhydrous	97.67
Potassium Chloride	400.00
Sodium Bicarbonate	3700.00
Sodium Chloride	6400.00
Sodium Phosphate, Monobasic, Monohydrate	125.00
<b>Specifications</b>	
Appearance	Clear, red solution
Aerobic bacteria	Not detected
Anaerobic bacteria	Not detected
Fungi (Yeast & Mold)	Not detected
pH at RT	6.9 - 7.5
Osmolality	260 - 340 mM/kg
Endotoxin	≤ 1.0 EU/ml
Mycoplasma	Not detected
Cell Culture - Cell Growth Promotion	3 Passages
Viability	≥ 75%

**THIS PRODUCT IS FOR LABORATORY USE ONLY.**