
	Document ID:	TDS-MCP-004-50L	Version:	006
	Date of Issue:	10-JAN-2025	Approved by:	Dr. Iman Kamranfar
	Review Date:	10-JAN-2027	Signature:	
	Title:	<b>TECHNICAL DATASHEET</b>		

### Dulbecco's Modified Eagle's Medium (DMEM)



<b>Filtration/ Treatment</b>	With high Glucose, L-Glutamine and without Sodium Pyruvate
<b>Product Code</b>	MCP-004-50L
<b>Shelf Life</b>	4 years from DOM
<b>Storage Temperature</b>	+2 to +8°C
<b>Shipping Temperature</b>	ambient

### QC Specifications

Physical and Chemical Analysis	Method	Specifications	Units
Appearance	Visual	Off-white to orange powder	n/a
Solubility	Internally validated	13.38 in ≤ 30 minutes	g/L
Moisture Content	Loss-on-drying	≤ 2.0 %	n/a
Particle Size (100 mesh)	Internally validated	≥ 90%	n/a
pH without NaHCO <sub>3</sub>	Electronic pH Meter	5.6 - 7.2	n/a
pH with NaHCO <sub>3</sub>	Electronic pH Meter	Record	n/a
Osmolality without NaHCO <sub>3</sub>	Osmometer	245 - 285	mOsm/kg
Osmolality with NaHCO <sub>3</sub>	Osmometer	Record	mOsm/kg
Endotoxin	LAL Kinetic	≤ 1.0	EU/ml
Mycoplasma	qPCR	Not detected	n/a
<b>Cell culture</b>			
Cell growth promotion	Passage test	3 passages	n/a
Viability	Trypan Blue Exclusion	≥75%	n/a

### Formulation

Amino Acids	CAS number	Concentration (mg/L)
Glycine	56-40-6	30.00
L-Arginine, Hydrochloride	1119-34-2	84.00
L-Cystine, Dihydrochloride	30925-07-6	63.00
L-Glutamine	56-85-9	584.00
L-Histidine, Hydrochloride, Monohydrate	5934-29-2	42.00
L-Isoleucine	73-32-5	105.00
L-Leucine	61-90-5	105.00
L-Lysine, Hydrochloride	657-27-2	146.00
L-Methionine	63-68-3	30.00
L-Phenylalanine	63-91-2	66.00
L-Serine	56-45-1	42.00
L-Threonine	72-19-5	95.00
L-Tryptophan	73-22-3	16.00
L-Tyrosine, Disodium, Dihydrate	122666-87-9	104.00
L-Valine	72-18-4	94.00
<b>Vitamins</b>		
Choline Chloride	67-48-1	4.00
D-Calcium Pantothenate	137-08-6	4.00
Folic Acid	59-30-3	4.00
Myo-Inositol	87-89-8	7.20
Niacinamide	98-92-0	4.00
Pyridoxine, Hydrochloride	58-56-0	4.00
Riboflavin	83-88-5	0.40
Thiamine, Hydrochloride	67-03-8	4.00
<b>Other Components</b>		
D-Glucose, Anhydrous	50-99-7	4500.00
Phenol Red, Sodium Salt	34487-61-1	15.00

	Document ID:	TDS-MCP-004-50L	Version:	006
	Date of Issue:	10-JAN-2025	Approved by:	Dr. Iman Kamranfar
	Review Date:	10-JAN-2027	Signature:	
	Title:	<b>TECHNICAL DATASHEET</b>		

Inorganic Salts		
Calcium Chloride, Anhydrous	10043-52-4	200.00
Ferric Nitrate, Nonahydrate	7782-61-8	0.10
Magnesium Sulfate, Anhydrous	7487-88-9	97.67
Potassium Chloride	7447-40-7	400.00
Sodium Chloride	7647-14-5	6400.00
Sodium Phosphate, Monobasic, Monohydrate	10049-21-5	125.00

Note: Sodium Bicarbonate is not included in the powdered media (3700.00 mg/L).

**Product Use:** This product is intended for laboratory use only.