

Title:	TECHNICAL DATASHEET		
Review Date:	10-JAN-2025	Signature:	Har
Date of Issue:	10-JAN-2023	Approved by:	Dr. Iman Kamranfar
Document ID:	TDS-S-FBS-SA-015	Version:	006

## **Foetal Bovine Serum**

Filtration, Treatment Triple 0.1µm sterile filtered		
Origin	Brazil	
Product Code	S-FBS-SA-015	
Shelf Life	5 years from DOM	
Storage Temperature	<-15°C	
Shipping Temperature	dry ice	

# **QC Specifications**

Physical and Chemical Analysis	Method	Specifications	Units
Appearance	Visual	Clear yellow-amber	n/a
Identity	Internally Validated	Bovine	n/a
pH at RT	Electronic pH Meter	6.8 - 8.2	n/a
Osmolality	Osmometer	260 - 340	mOsm/kg
Endotoxin	LAL Kinetic		EU/ml
Free Hemoglobin		< 10.0	
Specific Gravity	Colorimetric  Mass Balance	< 25.0	mg/dl g/ml
Specific Gravity	Mass Balance	> 1.01	g/mi
Sterility		I	
Sterility	Eur. Ph. 2.6.1	Pass	n/a
Mycoplasma	qPCR	Not detected	n/a
Wycopiasina	qrcn	Not detected	11/ a
Protein Profile			
Total Protein	IDEXX Catalyst One	3.0 - 4.5	g/dl
Albumin	IDEXX Catalyst One	1.4 - 3.4	g/dl
Globulin	IDEXX Catalyst One	0.4 - 2.4	g/dl
IgG	ELISA	< 400	μg/ml
Electrophoretic Pattern	Capillary Electrophoresis	Normal	n/a
zicoti opiioi otio i attein	capmany Electrophic colo		1.1,5
Antibiotic Testing			
Tetracycline	IDEXX Snap Test	Test and report	n/a
Oxytetracycline	IDEXX Snap Test	Test and report	n/a
Chlortetracycline	IDEXX Snap Test	Test and report	n/a
			1,72
Virus Testing	-	1	l .
BVDV/BHV-1/PI-3 (CPE)	Cell Culture	Not detected	n/a
Rabies Virus	qPCR	Not detected	n/a
Bluetongue Virus (BTV)	qPCR	Not detected	n/a
BRSV	qPCR	Not detected	n/a
Reo Virus	qPCR	Not detected	n/a
BAV	qPCR	Not detected	n/a
BoPV-1, -2	qPCR	Not detected	n/a
	•		
Antibody Testing			
	Serum Neutralization		
BVDV - Antibody Type 1	Test (Cell Culture) or Detection of Antibodies (ELISA)	Test and report	n/a
BVDV - Antibody Type 2	Serum Neutralization Test (Cell Culture) or Detection of Antibodies (ELISA)	Test and report	n/a
BHV-1	Detection of Antibodies (ELISA)	Test and report	n/a
PI-3	Detection of Antibodies (ELISA)	Test and report	n/a
Stark and the	T	T	T
Biochemistry	IDENA CALLAS		,
Aspartate Aminotransferase (AST)	IDEXX Catalyst One	Record	U/L
Alanine Aminotransferase (ALT)	IDEXX Catalyst One	Record	U/L
Lactate Dehydrogenase (LDH)	IDEXX Catalyst One	Record	U/L
Alkaline Phosphatase (ALKP)	IDEXX Catalyst One	Record	U/L
Gamma-Glutamyl Trans.(GGT)	IDEXX Catalyst One	Record	U/L



Document ID:  Date of Issue:	TDS-S-FBS-SA-015 10-IAN-2023	Version: Approved by:	006 Dr. Iman Kamranfar
Review Date:		Signature:	M
Title.			
Title:	TECHNICAL DATASHEET		

Cholesterol (CHOL)	IDEXX Catalyst One	Record	mmol/L
Glucose (GLU)	IDEXX Catalyst One	Record	mmol/L
Urea (BUN)	IDEXX Catalyst One	Record	mmol/L
Creatinine (CREA)	IDEXX Catalyst One	Record	μmol/L
Uric Acid (URIC)	IDEXX Catalyst One	Record	μmol/L
Calcium (Ca)	IDEXX Catalyst One	Record	mmol/L
Phosphorus (PHOS)	IDEXX Catalyst One	Record	mmol/L
Total Bilirubin (TBIL)	IDEXX Catalyst One	Record	μmol/L
Magnesium (Mg)	IDEXX Catalyst One	Record	mmol/L
Sodium (Na)	IDEXX Catalyst One	Record	mmol/L
Potassium (K)	IDEXX Catalyst One	Record	mmol/L
Chloride (Cl)	IDEXX Catalyst One	Record	mmol/L
Chioride (Ci)	IDEAX Catalyst Offe	Record	IIIIIIOI/L
Call Cultura Tastina Cution 1			
Cell Culture Testing - Option 1 Cell Line	Mathad	Succifications	Dogulto
Cell Line	Method	Specifications	Results
L-929, HeLa, MRC-5	Morphology	Tested vs. Control Serum	Scoring
			System 1
L-929, HeLa, MRC-5	Density	Tested vs. Control Serum	Scoring System 2
		Cell count [log10/ml]/dead cells vs.	
L-929, HeLa, MRC-5	Cell Count	Control	Record
Scoring system	Meaning		Results
	Dead Cells		0
	Many Cells degenerate		+ -
	and many dead cells		1
	Cells partially		
1 - Morphology	degenerate and many		2
	dead cells		1
	Few cells degenerate		
	and few dead cells		3
	Without pathological		
	findings		4
	Single cells/cell		
	aggregates		0
	< 50% confluency		1
2 - Density	50 - 90% confluency		2
	·		
	confluency		3
	overly confluent		4
Cell Culture Testing - Option 2 Cell Line	I advative d	Constitution	T
Cell Line	Method	Specifications	Units
	Multiple Passage -	>75% of control growth	
BHK-21, MRC-5	Records results		%
,	vs. control at day: 0, 3,		
	6, 12		1
	Plating Efficiency -	>75% PE vs. control PE	
BHK-21, MRC-5	Records results		%
	vs. control at day: 0, 3,		
	6, 12		
	Cloning Efficiency -		
BHK-21, MRC-5	Records results	>75% CE vs. control CE	%
	vs. control at day: 0, 3,		'-
	6, 12		



Document ID:	TDS-S-FBS-SA-015	Version:	006
Date of Issue:	10-JAN-2023	Approved by:	Dr. Iman Kamranfar
Review Date:	10-JAN-2025	Signature:	An
Title:	TECHNICAL DATASHEET		

#### **RECOMMENDED USE:**

#### Storage

To effectively preserve the integrity of animal serum, it should be stored frozen and protected from light. The recommended storage temperature is <-15°C.

Multiple thaw/freeze cycles should be avoided, as they will accelerate the degradation of serum nutrients and can encourage the formation of insoluble precipitates. For this reason, serum should never be stored in "frost-free" freezers. These types of freezers periodically warm themselves to avoid internal ice deposits and are detrimental to the stability of frozen serum products.

### **Suggested Thawing Procedure**

- 1. Remove the serum bottles from the freezer and allow them to adjust to room temperature for approximately 10 minutes.
- 2. Place each container in a 30 to 37 °C water bath or incubator. Excessive temperatures will degrade heat labile nutrients. If using a water bath, prevent the bottle caps from being submerged.
- 3. Gently agitate the bottles every 10 15 minutes until the serum is completely thawed.

#### **Efficient and Effective Usage**

After thawing, use the serum promptly. Liquid serum may be stored refrigerated (2 to 8 °C) up to four weeks. To avoid thaw/freeze cycles or long periods of refrigeration, it is recommended that any unused serum be immediately dispensed into small, useful aliquots and refrozen for future use.

Product use: This product is not intended for human or animal consumption or therapeutic use.